



More information on the website
radwag.com/en/info,w1,WRY

F1 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 5 kg), Plastic Box



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F1 |
| Construction | |
| Material | Stainless steel |
| Weight | 11 110 g |
| Nominal value | 1 g – 5 kg |
| Quantity | 16 pcs |

Compatible with

Tweezers

Brushes