




More information on the website
radwag.com/de/info,w1,X9P

Prüfgewicht der F2 Klasse - Schlitzgewicht- 200 g

OS-101-0028



The drawings, photos and graphics used are for illustrative purposes only.

Technische Daten

Messtechnische Parameter

| | |
|-------------|-------|
| OIML-Klasse | F2 |
| Nominalwert | 200 g |

Kompatibel mit (Additional Fee)

| | |
|-----------|--------|
| Pinzetten | Pinset |
|-----------|--------|