



# Knob weights, Set (1 g - 100 g), Plastic Box

OK-502-0005

More information on the website  
[radwag.com/en/info,w1,XHB](http://radwag.com/en/info,w1,XHB)



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

| Metrological parameters |                 |
|-------------------------|-----------------|
| OIML Class              | E1              |
| Construction            |                 |
| Material                | Stainless steel |
| Weight                  | 210 g           |
| Nominal value           | 1 – 100 g       |
| Quantity                | 9 pcs           |

## Compatible with (Additional Fee)

Tweezers

Brushes