



More information on the website  
[radwag.com/us/info,w1,XRY](http://radwag.com/us/info,w1,XRY)

# KR-01 Converter



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

### Physical parameters

Packaging dimensions	220×190×90 mm
Net weight	0,4 kg
Gross weight	0,8 kg

## Compatible with

MA X7.IC.A Moisture Analyzer  
 MA 5Y Moisture Analyzer  
 Tunnel detectors  
 C315 Platform scale  
 MA X2.IC.A Moisture Analyzers  
 AS X7 Analytical Balance  
 PS X2 Precision Balance  
 AS X2 PLUS Analytical Balances  
 C315.D Baby scale  
 WLC Precision Balance  
 C315.4P2 Beam Scales  
 WLC X7 Precision balance  
 PS X2 Balances for PGC  
 MA X2.A Moisture Analyzers

C315.0 Personal scale  
 PS X7 Precision Balance  
 C315.4P pallet scales  
 MA 5Y.IC Moisture Analyzer  
 AS 5Y Analytical Balance  
 H315 4C Platform Scale  
 C315.0W Personal scale  
 WLC C/2 Precision balance  
 PS 5Y Precision Balance  
 Postal scales for packages  
 PUE C315 Terminal  
 C315.4B Bed scale  
 PUE H315 Terminal  
 C315.K Chair scale

## Device dimensions

