



More information on the website
radwag.com/us/info,w1,Y0T

P0108.5 Cable

WX-013-0161



The drawings, photos and graphics used are for illustrative purposes only.

Overview

RS 232 cable for connecting a balance/scale to a computer.

Cable length: 5 m

Compatible with:

- Scales operated using the PUE C315 terminal
- PUE C315 terminal

Compatible with (Additional Fee)

MA X7.IC.A Moisture Analyzer
MUYA Modules
WPY Multifunctional Scale
C32 Multifunctional Scale
MA 5Y Moisture Analyzer
C315 Platform scale
WLY Scales for PGC
MA X2.IC.A Moisture Analyzer
AS X7 Analytical Balance
PS X2 Precision Balance
AS X2 PLUS Analytical Balance

PS 5Y Balances for PGC
AKM-2 Automatic Mass Comparator
PS R2.H Precision Balance
WTC Precision Balance
H315 4C Platform Scale
C315.OW Personal scale
Analytical Balance for Filters XA 5Y.F
WLY Precision Balance
UMA Automatic Mass Comparator
APP 5Y.KO Manual Mass Comparator
Robotic Mass Comparators

PS X7 Balances for PGC
C315.D Baby scale
WLC Precision Balance
PS R2 Precision Balance
C315.4P2 Beam Scales
WLC X7 Precision balance
PS X2 Balances for PGC
MA X2.A Moisture Analyzer
MAS Modules
MPS Modules
C315.4 Platform Scale
C32.PM Precision Balance
MA X7.A Moisture Analyzer
C315.O Personal scale
PS X7 Precision Balance
WPY Scales for PGC
C315.4P pallet scales
WAY 5Y.KO Manual Mass Comparator
MA 5Y.IC Moisture Analyzer
MA R Moisture Analyzer
AS 5Y Analytical Balance
AK-4 Automatic Mass Comparator

Density measurement
WLC C/2 Precision balance
PS R1 Precision Balance
MA X2.A Moisture Analyzers
5Y.PM.KB Manual Mass Comparator
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PS 5Y Precision Balance
Postal scales for packages
Automatic Vacuum Mass Comparators
MAS1 Modules
PMV 5Y Moisture Analyzer
XA 5Y.KO Manual Mass Comparator
PUE C315 Terminal
UYA 5Y.KO Manual Mass Comparator
AS R1 PLUS Analytical Balance
C315.4B Bed scale
Nano-Comparator NANO.AK-4.500.5Y
C315.K Chair scale
C315.8B Ramps bed scale