



More information on the website  
[radwag.com/it/info,w1,Y5R](http://radwag.com/it/info,w1,Y5R)

# Peso standard E2 - filo di acciaio inossidabile - 200 mg

OP-103-0060



The drawings, photos and graphics used are for illustrative purposes only.

## Specifiche

| Parametri metrologici |                 |
|-----------------------|-----------------|
| Classe OIML           | E2              |
| Costruzione           |                 |
| Materiale             | Stainless steel |
| Valore nominale       | 200 mg          |

## Compatibile con (Additional Fee)

Pinzette

Pennelli