



More information on the website
radwag.com/en/info,w1,Y5W

105 mm carbon fiber tweezers for mass standards up to 50 g



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Construction | |
|---------------------|--------------|
| Material | carbon fiber |
| Physical parameters | |
| Length | 105 mm |

Compatible with

- E1 Mass Standard - knob weights, sets (1 mg - 10 kg)
- F1 Mass Standard - stainless steel wires
- F1 Mass Standard - knob weights with adjustment chamber, sets
- F2 Mass Standard - stainless steel wires
- F2 Mass Standard - Carrier weights
- E2 Mass Standard - stainless steel wires
- E2 Mass Standard - knob weights, sets (1 mg - 10 kg)
- F2 Mass Standard - knob weights with adjustment chamber, sets
- E2 Mass Standard - cylindrical weights, sets (1 g - 10 kg)
- F2 Mass Standard - Slotted weights
- M1 Mass Standard - flat polygonal sheets
- F2 Mass Standard - flat polygonal sheets
- F1 Mass Standard - knob weights without adjustment chamber, set
- F1 Mass Standard - flat polygonal sheets

- F1 Mass Standard - knob weights with adjustment chamber
- E2 Mass Standard - flat polygonal sheets
- E2 Mass Standard - knob weights
- F1 Mass Standard - flat polygonal sheets, sets
- F1 Mass Standard - knob weights without adjustment chamber
- E1 Mass Standard - stainless steel wires
- F2 Mass Standard - knob weights with adjustment chamber
- F2 Mass Standard - flat polygonal sheets, sets
- M1 Mass Standard - knob weights
- E1 Mass Standard - flat polygonal sheets
- E1 Mass Standard - knob weights
- E2 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
- E1 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
- F2 Mass Standard - Stainless steel wires, sets

E1 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)
F1 Mass Standard - cylindrical weights
E2 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)

F1 Mass Standard - stainless steel wires, sets
M1 Mass Standard - flat polygonal sheets, sets