



More information on the website
radwag.com/en/info,w1,YAN

R-Barcode

WX-010-0109



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Components and software

Operating system

Windows 7, 8, 10

Compatible with (Additional Fee)

MUYA Modules

HX7 H Stainless steel Multifunctional Scale

MAS1 Modules

HX7.4N H Multifunctional Stainless Steel Ramp Scale

HX7.4C Multifunctional Platform Scale

HX7 HR Multifunctional Scale

XA 5Y.M.A.P Microbalance for pipettes calibration

WPY Multifunctional Scale

MYA 5Y.P Microbalance for pipettes calibration

C32 Multifunctional Scale

HX5.EX-1.4P2 H Stainless Steel Beam Scales

HY10.HRP.H Stainless steel High Resolution Scale

MYA 5Y.F Microbalance for Filters

HX5.EX-1 HR Stainless Steel 1-Load-Cell Platform Scales

HX7.4P2 H Multifunctional Stainless Steel Beam Scales

WLY Scales for PGC

AS X7 Analytical Balance

PS X7 Balances for PGC

HX7 Multifunctional Scale

HX5.EX-1.4 H/Z Stainless Steel Platform Scales, Pit Version

HX7.4 H/Z Multifunctional Stainless Steel Platform Scale, pit version

PS X7 Precision Balance

WPY Scales for PGC

XA 5Y.A Analytical Balance

HX7.4 H Multifunctional Stainless Steel Platform Scale

HX7.4P2 Multifunctional Beam Scales

HX5.EX-1 C 1-Load-Cell Platform Scales

HX5.EX-1 H Waterproof Platform Scales

AS 5Y Analytical Balance

PS 5Y Balances for PGC

PS R2.H Precision Balance

PUE HX5.EX Weighing Terminal

HY10.H Scales for PGC

HX5.EX-1.4 H Stainless Steel 4-Load-Cell Platform Scales

CY10 Precision Balance

XA 5Y.F Analytical Balance for Filters

MWMH Modules
XA 5Y.M.A Microbalance
PS R2 Precision Balance
HX7.4N H.LD Multifunctional Stainless Steel Ramp Scale
HX5.EX-1.4P2 C Beam Scales
PUE CY10 Weighing Terminal
CY10 Multifunctional Scale
HX7.4P H Multifunctional Stainless Steel Pallet Scales
HX5.EX-1.4N Stainless Steel Ramp Scales
HY10 Multifunctional Scale
HX5.EX-1.4 C 4-Load-Cell Platform Scales
MAS Modules
HY10 H Stainless Steel Multifunctional Scale
MPS Modules
HX5.EX-1.4P H Stainless Steel Pallet Scales
HY10.HRP High Resolution Scale
C32.PM Precision Balance
PUE 7.1.HRP High Resolution Scale
Automatic Balances for. Pipettes Calibration

WLY Precision Balance
MW-01 Mass Converters
MYA 5Y Microbalance
PS R1 Precision Balance
XA 5Y Analytical Balance
AS R2 PLUS Analytical Balance
PUE C32 Weighing Terminal
PUE 7.1 Weighing Terminal
PUE HX7 Weighing Terminal
Barcode scanners
PS 5Y Precision Balance
HX5.EX-1.4P C Pallet Scales
XA 5Y.M Microbalance
PUE HY10 Weighing Terminal
AS R1 PLUS Analytical Balance
H315 H Waterproof Platform Scale
HX7.4P Multifunctional pallet scales
XA 5Y.M.A.S Microbalance for stents weighing