



More information on the website  
[radwag.com/en/info,w1,YG8](http://radwag.com/en/info,w1,YG8)

# SA/APP/H Antivibration Table for Mass Comparators stainless steel

WX-002-0040



The drawings, photos and graphics used are for illustrative purposes only.

## Datasheet

Construction	
Construction	Stainless steel AISI304
Load-bearing capacity	150 kg
Physical parameters	
Packaging dimensions	710×470×225 mm
Net weight	128 kg
Gross weight	142 kg
Length	710 mm
Width	470 mm
Height	225 mm

## Compatible with (Additional Fee)

APP 5Y.KO Manual Mass Comparator

5Y.PM.KB Manual Mass Comparator

## Device dimensions

