



More information on the website
radwag.com/en/info,w1,YJD

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 g - 2 kg), Plastic Box

OK-302-0029



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F2 |
| Construction | |
| Material | Stainless steel |
| Weight | 6 110 g |
| Nominal value | 1 g – 2 kg |
| Quantity | 15 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes