



More information on the website
radwag.com/en/info,w1,YRW

SAL/C/PLUS Antivibration Table for PLUS Series Balances



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

Construction	
Construction	mild steel St3S
Load-bearing capacity	40 kg
Physical parameters	
Packaging dimensions	820×711×816 mm
Net weight	68 kg (table 28 kg + stone 40 kg)
Gross weight	77.5 kg
Length	770 mm
Width	680 mm
Height	815±10 mm

Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 Microbalance for Pipette Calibration MYA 5Y.P
 Microbalance for Filters MYA 5Y.F
 XA 5Y.M.A Microbalance

XA 5Y.F Analytical Balance for Filters
 MYA 5Y Microbalance
 XA 5Y Analytical Balance
 XA 5Y.M Microbalance

Device dimensions

SAL/C/PLUS

