



More information on the website
radwag.com/us/info,w1,Z4V

Under-Pan Weighing Hook (XA, AS, PS)



The drawings, photos and graphics used are for illustrative purposes only.

Overview

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

Compatible with

XA 5Y.M.A.P Microbalance for pipettes calibration
 AS X7 Analytical Balance
 PS X2 Precision Balance
 AS X2 PLUS Analytical Balances
 PS X7 Balances for PGC
 XA 5Y.M.A Microbalance
 PS R2 Precision Balance
 PS X2 Balances for PGC
 PS X7 Precision Balance
 XA 5Y.A Analytical Balance
 AS 5Y Analytical Balance

PS 5Y Balances for PGC
 PS R2.H Precision Balance
 Analytical Balance for Filters XA 5Y.F
 PS R1 Precision Balance
 XA 5Y Analytical Balance
 AS R2 PLUS Analytical Balance
 PS 5Y Precision Balance
 XA 5Y.M Microbalance
 AS R1 PLUS Analytical Balance
 XA 5Y.M.A.S Microbalance for stents weighing