



More information on the website  
[radwag.com/us/info,w1,Z4V](http://radwag.com/us/info,w1,Z4V)

# Under-Pan Weighing Hook (XA, AS, PS)

WX-015-0325



The drawings, photos and graphics used are for illustrative purposes only.

## Overview

A hook that allows under-pan weighing. It is compatible with XA, AS and PS series balances.

## Compatible with (Additional Fee)

XA 5Y.M.A.P Microbalance for pipettes calibration  
 AS X7 Analytical Balance  
 PS X2 Precision Balance  
 AS X2 PLUS Analytical Balances  
 PS X7 Balances for PGC  
 XA 5Y.M.A Microbalance  
 PS R2 Precision Balance  
 PS X2 Balances for PGC  
 PS X7 Precision Balance  
 XA 5Y.A Analytical Balance  
 AS 5Y Analytical Balance

PS 5Y Balances for PGC  
 PS R2.H Precision Balance  
 Analytical Balance for Filters XA 5Y.F  
 PS R1 Precision Balance  
 XA 5Y Analytical Balance  
 AS R2 PLUS Analytical Balance  
 PS 5Y Precision Balance  
 XA 5Y.M Microbalance  
 AS R1 PLUS Analytical Balance  
 XA 5Y.M.A.S Microbalance for stents weighing