



More information on the website
radwag.com/en/info,w1,ZF0

F2 Mass Standard - Knob Weights With Adjustment Chamber, Set (1 mg - 1 kg), Wooden Box

OK-301-0041



The drawings, photos and graphics used are for illustrative purposes only.

Datasheet

| Metrological parameters | |
|-------------------------|-----------------|
| OIML Class | F2 |
| Construction | |
| Material | Stainless steel |
| Weight | 2 111.11 g |
| Nominal value | 1 mg – 1 kg |
| Quantity | 25 pcs |

Compatible with (Additional Fee)

Tweezers

Brushes