



More information on the website  
[radwag.com/en/info,w1,ZHC](http://radwag.com/en/info,w1,ZHC)

## 200 mm tweezers with carbon fiber tip for mass standards up to 200 g



The drawings, photos and graphics used are for illustrative purposes only.

### Datasheet

Construction	
Material	Stainless steel, carbon fiber
Physical parameters	
Length	200 mm

### Compatible with

- E1 Mass Standard - knob weights, sets (1 mg - 10 kg)
- F1 Mass Standard - stainless steel wires
- F1 Mass Standard - knob weights with adjustment chamber, sets
- F2 Mass Standard - stainless steel wires
- F2 Mass Standard - Carrier weights
- E2 Mass Standard - stainless steel wires
- E2 Mass Standard - knob weights, sets (1 mg - 10 kg)
- F2 Mass Standard - knob weights with adjustment chamber, sets
- E2 Mass Standard - cylindrical weights, sets (1 g - 10 kg)
- F2 Mass Standard - Slotted weights
- M1 Mass Standard - flat polygonal sheets
- F2 Mass Standard - flat polygonal sheets
- F1 Mass Standard - knob weights without adjustment chamber, set
- F1 Mass Standard - flat polygonal sheets
- E2 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)
- F1 Mass Standard - knob weights with adjustment chamber
- E2 Mass Standard - knob weights
- E2 Mass Standard - flat polygonal sheets
- F1 Mass Standard - flat polygonal sheets, sets
- E1 Mass Standard - stainless steel wires
- F2 Mass Standard - knob weights with adjustment chamber
- F2 Mass Standard - flat polygonal sheets, sets
- M1 Mass Standard - knob weights
- E1 Mass Standard - flat polygonal sheets
- E1 Mass Standard - knob weights
- E2 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
- E1 Mass Standard - flat polygonal sheets, sets (1 mg - 500 mg)
- F2 Mass Standard - Stainless steel wires, sets

E1 Mass Standard - stainless steel wires, sets (1 mg - 500 mg)  
F1 Mass Standard - cylindrical weights  
F1 Mass Standard - knob weights without adjustment chamber

F1 Mass Standard - stainless steel wires, sets  
M1 Mass Standard - flat polygonal sheets, sets